

These vectors are provided with a glycerol stock of bacterial strain JM109.

Description: The pGEM®-3Z(+/-) Vectors[®] are derived from the pGEM®-3Z Vector and contain the origin of replication of the filamentous phage f1. These plasmids serve as standard cloning vectors, as templates for in vitro transcription, and can be used for the production of circular ssDNA.

The pGEM®-3Z(+) and pGEM®-3Z(-) Vectors are identical except for the orientation of the f1 origin.

Feature(s)

- **Blue/White Screening:** Allows the easy identification of recombinant clones.
- **Versatile:** These vectors can be used for standard cloning, single-stranded DNA production and in vitro transcription from SP6 and T7 RNA polymerase promoters flanking the multiple cloning region.
- **Convenient:** Multiple cloning site provides a selection of restriction sites for cloning.

Storage Conditions: Store vector at -20°C or -70°C and host strain at -70°C.

GenBank®/EMBL Accession Number: (+) X65306; (-) X65307.


Related Products

	Pg.
Riboprobe® Transcription Systems	5.3
Wizard® DNA Purification Systems	2.4
JM109 Competent Cells	13.11

Additional Information Available

	Lit.#
pGEM®-3Z(+) Sequence & Map	TB086
pGEM®-3Z(-) Sequence & Map	TB045

©U.S. Pat. No. 4,766,072.



For detailed vector maps, please visit our web site at www.promega.com

To: Bao-Qun Li
Tel: 703-305-1695
Fax: 703-746-7454

These vectors are provided with a glycerol stock of JM109.

Description: The pGEM®-5Z(+) and pGEM®-5Z(-) Vectors[®] are derived from the pGEM®-5Z Vector and contain the origin of replication f1. These plasmids serve as standard cloning vectors, as templates for in vitro transcription, and can be used for the production of circular ssDNA. The multiple cloning region contains T7 and SP6 RNA polymerase promoters flanking the α-peptide coding region of β-galactosidase (the α-peptide allows recombinant clones to be directly on indicator plates). The multiple cloning region contains: Apa I, Aat II, Sph I, Nco I, Sac II, EcoR V, Spe I, Not I, Pst I, and Nsi I. This arrangement is designed specifically for use with Promega's Erase-a-Base® System. The pGEM Vectors are identical except for the orientation of the f1 origin.

Feature(s)

- **Blue/White Screening:** Allows the easy identification of recombinant clones.
- **Versatile:** These vectors can be used for standard cloning, single-stranded DNA production and in vitro transcription from SP6 and T7 RNA polymerase promoters flanking the multiple cloning region.
- **Convenient:** Multiple cloning site provides a selection of restriction sites for cloning.
- **Unidirectional Deletions:** Restriction sites for use with Promega's Erase-a-Base® System

Storage Conditions: Store vector at -20°C or -70°C.

GenBank®/EMBL Accession Number: (+) X65306; (-) X65307.

Reference(s)

1. Yanisch-Perron, C., Vieira, J. and Messing, J.


Related Products

Riboprobe® Transcription Systems
Wizard® DNA Purification Systems
Erase-a-Base® System
pGEM®-T Vectors Systems
JM109 Competent Cells

Additional Information Available

pGEM®-5Z(+) Sequence & Map
pGEM®-5Z(-) Sequence & Map

©U.S. Pat. No. 4,766,072.



For detailed vector maps, please visit our web site at www.promega.com

TB 086 (+)
TB 045 (-)

The pGEM®-3Z Vector is provided with a glycerol stock of bacterial strain JM109.

Description: The pGEM®-3Z Vector^(a) is intended for use as a standard cloning vector, as well as for the highly efficient synthesis of RNA in vitro. The vector carries the *lacZ* α-peptide and the multiple cloning region arrangement from pUC18 (1). In addition, the vector contains both the SP6 and T7 RNA polymerase promoters flanking the multiple cloning region.

The pGEM®-3Z and pGEM®-4Z Vectors are identical except for the orientation of the SP6 and T7 promoters.

Feature(s)

- **Blue/White Screening:** Allows the easy identification of recombinant clones.
- **Versatile:** This vector can be used for standard cloning and in vitro transcription from SP6 and T7 RNA polymerase promoters flanking the multiple cloning region.
- **Convenient:** Multiple cloning site provides a selection of restriction sites for cloning.

Storage Conditions: Store vector at -20°C or -70°C and host strain at -70°C.

GenBank®/EMBL Accession Number: X65304.

Reference(s)

1. Yanisch-Perron, C. *et al.* (1985) *Gene* **33**, 103.


Related Products

	Pg.
Riboprobe® Transcription Systems	5.3
Wizard® DNA Purification Systems	2.4
JM109 Competent Cells	13.11

Additional Information Available

	Lit. #
Sequence & Map	TB033

©U.S. Pat. No. 4,766,072.



The pGEM®-4Z Vector is provided with a glycerol stock of bacterial strain JM109.

Description: The pGEM®-4Z Vector^(a) is inter cloning vector, as well as for the highly efficient vector carries the *lacZ* α-peptide and the multi from pUC18 (1) allowing recombinants to be screened. In addition, the vector contains both polymerase promoters flanking the multiple cloning region.

The pGEM®-3Z and pGEM®-4Z Vectors are identical except for the orientation of the SP6 and T7 promoters.

Feature(s)

- **Blue/White Screening:** Allows the easy identification of recombinant clones.
- **Versatile:** This vector can be used for standard cloning and in vitro transcription from SP6 and T7 RNA polymerase promoters flanking the multiple cloning region.
- **Convenient:** Multiple cloning site provides a selection of restriction sites for cloning.

Storage Conditions: Store vector at -20°C or -70°C and host strain at -70°C.

GenBank®/EMBL Accession Number: X65304.

Reference(s)

1. Yanisch-Perron, C. *et al.* (1985) *Gene* **33**, 103.


Related Products

	Pg.
Riboprobe® Transcription Systems	5.3
Wizard® DNA Purification Systems	2.4
JM109 Competent Cells	13.11

Additional Information Available

	Lit. #
Sequence & Map	TB033


©U.S. Pat. No. 4,766,072.



For detailed vector maps, please visit our website.

P2271
pGEM


Close Figure Window
Print Page



Web Figure Number: 0282vc
Figure Display Window

P2271
pGEM

Close Figure Window
Print Page



Web Figure Number: 0282vc
Figure Display Window

pGEM®*-3Zf(+) vector sequence

This vector can be obtained from Promega Corporation, Madison WI.
Call one of the following numbers for order or technical information.

Order or Technical 1-800-356-9526
Outside U.S. 608-274-4330

pGEM®-3Zf(+) vector sequence reference points:

Base pairs	3197
T7 RNA transcription initiation site	1
SP6 RNA transcription initiation site	69
T7 RNA polymerase promoter (-17 to +3)	3181-3
SP6 RNA polymerase promoter (-17 to +3)	67-86
multiple cloning region	5-61
phage f1 region	2562-3017
lacZ start codon	108
lac operon sequences	3018-3178; 94-323
lac operator	128-144
beta-lactamase (Amp ^r) coding region	1265-2125
binding site of pUC/M13 Forward Sequencing Primer	3138-3154
binding site of pUC/M13 Reverse Sequencing Primer	104-120

pGEM is a registered trademark of Promega Corporation, Madison, WI.

*U.S. Pat. No. 4,766,072 has been issued to Promega Corporation for transcription ve
different bacteriophage RNA polymerase promoter sequences separated by a series of u
sites into which foreign DNA can be inserted.

```
1   GGGCGAATTC GAGCTCGGTA CCCGGGGATC CTCTAGAGTC GACCTGCAGG
51  CATGCAAGCT TGAGTATTCT ATAGTGTAC CTAATAGCT TGGCGTAATC
101 ATGGTCATAG CTGTTTCCTG TGTGAAATG TTATCCGCTC ACAATTCCAC
151 ACAACATACG AGCCGGAAGC ATAAAGTGTA AAGCCTGGGG TGCCTAATGA
201 GTGAGCTAAC TCACATTAAT TCGGTTGCGC TCACTGCCCC CTTTCCAGTC
251 GGGAAACCTG TCGTGCCAGC TGCATTAATG AATCGGCCAA CGCGCGGGGA
301 GAGGCGGTTT GCGTATTGGG CGCTCTTCCG CTTCTCGCT CACTGACTCG
351 CTGCGCTCGG TCGTTCGGCT GCGGCGAGCG GTATCAGCTC ACTCAAAGGC
401 GGTAAATACG TTATCCACAG AATCAGGGGA TAACGCAGGA AAGAACATGT
451 GAGCAAAAGG CCAGCAAAAG GCCAGGAACC GTAAAAAGGC CGCGTTGCTG
501 GCGTTTTTCC ATAGGCTCCG CCCCCCTGAC GAGCATCACA AAAATCGACG
551 CTCAAGTCAG AGGTGGCGAA ACCCGACAGG ACTATAAAGA TACCAGGCGT
601 TTCCCCCTGG AAGCTCCCTC GTGCGCTCTC CTGTTCCGAC CCTGCCGCTT
651 ACCGGATACC TGTCCGCCTT TCTCCCTTCG GGAAGCGTGG CGCTTTCTCA
701 TAGCTCACGC TGTAGGTATC TCAGTTCGGT GTAGGTCGTT CGCTCCAAGC
751 TGGGCTGTGT GCACGAACCC CCCGTTACGC CCGACCGCTG CGCCTTATCC
801 GGTAACATATC GTCTTGAGTC CAACCCGGTA AGACACGACT TATCGCCACT
```

851 GGCAGCAGCC ACTGGTAACA GGATTAGCAG AGCGAGGTAT GTAGGCGGTG
901 CTACAGAGTT CTTGAAGTGG TGGCCTAACT ACGGCTACAC TAGAAGAACA
951 GTATTTGGTA TCTGCGCTCT GCTGAAGCCA GTTACCTTCG GAAAAAGAGT
1001 TGGTAGCTCT TGATCCGGCA AACAAACCAC CGCTGGTAGC GGTGGTTTTT
1051 TTGTTTGCAA GCAGCAGATT ACGCGCAGAA AAAAAGGATC TCAAGAAGAT
1101 CCTTTGATCT TTTCTACGGG GTCTGACGCT CAGTGAACG AAAACTCACG
1151 TTAAGGGATT TTGGTCATGA GATTATCAAA AAGGATCTTC ACCTAGATCC
1201 TTTTAAATTA AAAATGAAGT TTTAAATCAA TCTAAAGTAT ATATGAGTAA
1251 ACTTGGTCTG ACAGTTACCA ATGCTTAATC AGTGAGGCAC CTATCTCAGC
1301 GATCTGTCTA TTTCGTTTAT CCATAGTTGC CTGACTCCCC GTCGTGTAGA
1351 TAACTACGAT ACGGGAGGGC TTACCATCTG GCCCCAGTGC TGCAATGATA
1401 CCGCGAGACC CACGCTCACC GGCTCCAGAT TTATCAGCAA TAAACCAGCC
1451 AGCCGGAAGG GCCGAGCGCA GAAGTGGTCC TGCAACTTTA TCCGCCTCCA
1501 TCCAGTCTAT TAATTGTTGC CGGGAAGCTA GAGTAAGTAG TTCGCCAGTT
1551 AATAGTTTGC GCAACGTTGT TGCCATTGCT ACAGGCATCG TGGTGTACG
1601 CTCGTCGTTT GGTATGGCTT CATTGAGCTC CGGTTCCCAA CGATCAAGGC
1651 GAGTTACATG ATCCCCCATG TTGTGCAAAA AAGCGGTTAG CTCCTTCGGT
1701 CCTCCGATCG TTGTCAGAAG TAAGTTGGCC GCAGTGTAT CACTCATGGT
1751 TATGGCAGCA CTGCATAATT CTCTTACTGT CATGCCATCC GTAAGATGCT
1801 TTTCTGTGAC TGGTGAGTAC TCAACCAAGT CATTCTGAGA ATAGTGTATG
1851 CGGCGACCGA GTTGCTCTTG CCCGGCGTCA ATACGGGATA ATACCGCGCC
1901 ACATAGCAGA ACTTTAAAAG TGCTCATCAT TGGAAAACGT TCTTCGGGGC
1951 GAAAACTCTC AAGGATCTTA CCGCTGTTGA GATCCAGTTC GATGTAACCC
2001 ACTCGTGCAC CCAACTGATC TTCAGCATCT TTTACTTTCA CCAGCGTTTC
2051 TGGGTGAGCA AAAACAGGAA GGCAAAATGC CGCAAAAAG GGAATAAGGG
2101 CGACACGGAA ATGTTGAATA CTCATACTCT TCCTTTTTC AATTATTGA
2151 AGCATTTATC AGGGTTATTG TCTCATGAGC GGATACATAT TTGAATGTAT
2201 TTAGAAAAAT AAACAAATAG GGGTTCCGCG CACATTTCCC CGAAAAGTGC
2251 CACCTGACGT CTAAGAAACC ATTATTATCA TGACATTAAC CTATAAAAAT
2301 AGGCGTATCA CGAGGCCCTT TCGTCTCGCG CGTTTCGGTG ATGACGGTGA
2351 AAACCTCTGA CACATGCAGC TCCCGGAGAC GGTACAGCT TGTCTGTAAG
2401 CGGATGCCGG GAGCAGACAA GCCCGTCAGG GCGCGTCAGC GGGTGTGGC

2451 GGGTGTCTGGG GCTGGCTTAA CTATGCGGCA TCAGAGCAGA TTGTACTGAG
2501 AGTGCACCAT ATGCGGTGTG AAATACCGCA CAGATGCGTA AGGAGAAAAT
2551 ACCGCATCAG GAAATTGTAA GCGTTAATAT TTTGTTAAAA TTCGCGTTAA
2601 ATTTTGTGTA AATCAGCTCA TTTTTTAACC AATAGGCCGA AATCGGCAAA
2651 ATCCCTTATA AATCAAAAGA ATAGACCGAG ATAGGGTTGA GTGTTGTTCC
2701 AGTTTGAAC AAGAGTCCAC TATTAAAGAA CGTGGACTCC AACGTCAAAG
2751 GCGGAAAAAC CGTCTATCAG GGCGATGGCC CACTACGTGA ACCATCACCC
2801 TAATCAAGTT TTTTGGGGTC GAGGTGCCGT AAAGCACTAA ATCGGAACCC
2851 TAAAGGGAGC CCCCATTATA GAGCTTGACG GGGAAAGCCG GCGAACGTGG
2901 CGAGAAAGGA AGGGAAGAAA GCGAAAGGAG CGGGCGCTAG GCGCTGGCA
2951 AGTGTAGCGG TCACGCTGCG CGTAACCACC ACACCGCCG CGCTTAATGC
3001 GCCGCTACAG GGCGCGTCCA TTCGCCATTC AGGCTGCGCA ACTGTTGGGA
3051 AGGGCGATCG GTGCGGGCCT CTTGCTATT ACGCCAGCTG GCGAAAGGGG
3101 GATGTGCTGC AAGGCGATTA AGTTGGGTAA CGCCAGGGTT TTCCCAGTCA
3151 CGACGTTGTA AAACGACGGC CAGTGAATTG TAATACGACT CACTATA

Sequence and reference points updated 06-Jul-99.

pGEM®*-3Zf(-) Vector Sequence

This vector can be obtained from Promega Corporation, Madison WI.
Call one of the following numbers for order or technical information.

Order or Technical 1-800-356-9526
Outside U.S. 608-274-4330

pGEM®-3Zf(-) Vector sequence reference points:

Base pairs	3197
T7 RNA polymerase transcription initiation site	1
SP6 RNA polymerase transcription initiation site	69
T7 RNA polymerase promoter (-17 to +3)	3181-3
SP6 RNA polymerase promoter (-17 to +3)	67-86
Multiple cloning region	5-61
Phage f1 region	2562-3017
lacZ start codon	108
lac operon sequences	3018-3178; 94-323
beta-lactamase coding region	1265-2125
Binding site of pUC/M13 forward sequencing primer	3138-3154
Binding site of pUC/M13 reverse sequencing primer	104-120

pGEM is a registered trademark of Promega Corporation, Madison, WI.

*U.S. Pat. No. 4,766,072 has been issued to Promega Corporation for transcription vectors having two different bacteriophage RNA polymerase promoter sequences separated by a series of unique restriction sites into which foreign DNA can be inserted.

..

```
1  GGGCGAATTC GAGCTCGGTA CCCGGGGATC CTCTAGAGTC GACCTGCAGG
51  CATGCAAGCT TGAGTATTCT ATAGTGTAC CTAATAGCT TGGCGTAATC
101 ATGGTCATAG CTGTTTCCTG TGTGAAATG TTATCCGCTC ACAATTCAC
151 ACAACATACG AGCCGGAAGC ATAAAGTGT AAGCCTGGGG TGCCTAATGA
201 GTGAGCTAAC TCACATTAAT TCGTTGCGC TCACTGCCCC CTTTCCAGTC
251 GGGAAACCTG TCGTGCCAGC TGCATTAATG AATCGGCCAA CGCGCGGGGA
301 GAGGCGGTTT GCGTATTGGG CGCTCTTCCG CTTCTCGCT CACTGACTCG
351 CTGCGCTCGG TCGTTCGGCT GCGGCGAGCG GTATCAGCTC ACTCAAAGGC
401 GGTAAATACG TTATCCACAG AATCAGGGGA TAACGCAGGA AAGAACATGT
451 GAGCAAAAGG CCAGCAAAAG GCCAGGAACC GTAAAAAGGC CGCGTTGCTG
501 GCGTTTTTCC ATAGGCTCCG CCCCCCTGAC GAGCATCACA AAAATCGACG
551 CTCAAGTCAG AGGTGGCGAA ACCCGACAGG ACTATAAAGA TACCAGGCGT
601 TTCCCCCTGG AAGCTCCCTC GTGCGCTCTC CTGTTCCGAC CCTGCCGCTT
651 ACCGGATACC TGTCCGCCTT TCTCCCTTCG GGAAGCGTGG CGCTTTCTCA
701 TAGCTCACGC TGTAGGTATC TCAGTTCGGT GTAGGTCGTT CGCTCCAAGC
751 TGGGCTGTGT GCACGAACCC CCCGTTACAG CCGACCGCTG CGCCTTATCC
801 GGTAACATATC GTCTTGAGTC CAACCCGGTA AGACACGACT TATCGCCACT
```


851 GGCAGCAGCC ACTGGTAACA GGATTAGCAG AGCGAGGTAT GTAGGCGGTG
901 CTACAGAGTT CTTGAAGTGG TGGCCTAACT ACGGCTACAC TAGAAGAACA
951 GTATTTGGTA TCTGCGCTCT GCTGAAGCCA GTTACCTTCG GAAAAAGAGT
1001 TGGTAGCTCT TGATCCGGCA AACAAACCAC CGCTGGTAGC GGTGGTTTTT
1051 TTGTTTGCAA GCAGCAGATT ACGCGCAGAA AAAAAGGATC TCAAGAAGAT
1101 CCTTTGATCT TTTCTACGGG GTCTGACGCT CAGTGAACG AAAACTCACG
1151 TTAAGGGATT TTGGTCATGA GATTATCAAA AAGGATCTTC ACCTAGATCC
1201 TTTTAAATTA AAAATGAAGT TTTAAATCAA TCTAAAGTAT ATATGAGTAA
1251 ACTTGGTCTG ACAGTTACCA ATGCTTAATC AGTGAGGCAC CTATCTCAGC
1301 GATCTGTCTA TTTCGTTTCAT CCATAGTTGC CTGACTCCCC GTCGTGTAGA
1351 TAACTACGAT ACGGGAGGGC TTACCATCTG GCCCCAGTGC TGCAATGATA
1401 CCGCGAGACC CACGCTCACC GGCTCCAGAT TTATCAGCAA TAAACCAGCC
1451 AGCCGGAAGG GCCGAGCGCA GAAGTGGTCC TGCAACTTTA TCCGCCTCCA
1501 TCCAGTCTAT TAATTGTTGC CGGGAAGCTA GAGTAAGTAG TTCGCCAGTT
1551 AATAGTTTGC GCAACGTTGT TGCCATTGCT ACAGGCATCG TGGTGTACAG
1601 CTCGTCGTTT GGTATGGCTT CATTTCAGCTC CGGTTCCCAA CGATCAAGGC
1651 GAGTTACATG ATCCCCCATG TTGTGCAAAA AAGCGGTTAG CTCCTTCGGT
1701 CCTCCGATCG TTGTCAGAAG TAAGTTGGCC GCAGTGTAT CACTCATGGT
1751 TATGGCAGCA CTGCATAATT CTCTTACTGT CATGCCATCC GTAAGATGCT
1801 TTTCTGTGAC TGGTGAGTAC TCAACCAAGT CATTCTGAGA ATAGTGTATG
1851 CGGCGACCGA GTTGCTCTTG CCCGGCGTCA ATACGGGATA ATACCGCGCC
1901 ACATAGCAGA ACTTTAAAAG TGCTCATCAT TGGAAAACGT TCTTCGGGGC
1951 GAAAACTCTC AAGGATCTTA CCGCTGTTGA GATCCAGTTC GATGTAACCC
2001 ACTCGTGCAC CCAACTGATC TTCAGCATCT TTTACTTTCA CCAGCGTTTC
2051 TGGGTGAGCA AAAACAGGAA GGCAAAATGC CGCAAAAAG GGAATAAGGG
2101 CGACACGGAA ATGTTGAATA CTCATACTCT TCCTTTTTCA ATATTATTGA
2151 AGCATTTATC AGGGTTATTG TCTCATGAGC GGATACATAT TTGAATGTAT
2201 TTAGAAAAAT AAACAAATAG GGGTTCCGCG CACATTTCCC CGAAAAGTGC
2251 CACCTGACGT CTAAGAAACC ATTATTATCA TGACATTAAC CTATAAAAT
2301 AGGCGTATCA CGAGGCCCTT TCGTCTCGCG CGTTTCGGTG ATGACGGTGA
2351 AAACCTCTGA CACATGCAGC TCCCGGAGAC GGTCACAGCT TGTCTGTAAG
2401 CGGATGCCGG GAGCAGACAA GCCCCTCAGG GCGCGTCAGC GGGTGTGGC

2451 GGGTGTCCGG GCTGGCTTAA CTATGCGGCA TCAGAGCAGA TTGTACTGAG
2501 AGTGCACCAT ATGCGGTGTG AAATACCGCA CAGATGCGTA AGGAGAAAAT
2551 ACCGCATCAG GACGCGCCCT GTAGCGGCGC ATTAAGCGCG GCGGGTGTGG
2601 TGGTTACGCG CAGCGTGACC GCTACACTTG CCAGCGCCCT AGCGCCCGCT
2651 CCTTTCGCTT TCTTCCTTC CTTTCTCGCC ACGTTCGCCG GCTTTCCTCG
2701 TCAAGCTCTA AATCGGGGGC TCCCTTTAGG GTTCCGATTT AGTGCTTTAC
2751 GGCACCTCGA CCCCCAAAAA CTTGATTAGG GTGATGGTTC ACGTAGTGGG
2801 CCATCGCCCT GATAGACGGT TTTTCGCCCT TTGACGTTGG AGTCCACGTT
2851 CTTTAATAGT GGACTCTTGT TCCAAACTGG AACAACTC AACCTATCT
2901 CGGTCTATTC TTTTGATTTA TAAGGGATTT TGCCGATTTT GGCCTATTGG
2951 TTAAAAAATG AGCTGATTTA AAAAAATTT AACGCGAATT TTAACAAAAT
3001 ATTAACGCTT ACAATTTCCA TTCGCCATTC AGGCTGCGCA ACTGTTGGGA
3051 AGGGCGATCG GTGCGGGCCT CTTGCTATT ACGCCAGCTG GCGAAAGGGG
3101 GATGTGCTGC AAGGCGATTA AGTTGGGTAA CGCCAGGGTT TTCCAGTCA
3151 CGACGTTGTA AAACGACGGC CAGTGAATTG TAATACGACT CACTATA

Vector sequence updated 17-May-00.